



VaQua™ Coated Intermittent Catheter
Activates in record time

 **COVIDIEN**
positive results for life™

COVIDIEN, COVIDIEN with logo and ™ marked brands are trademarks of Covidien AG or its affiliate. © 2007 Covidien AG or its affiliate. All rights reserved. Covered by one or more of US Patent Nos. 6, 645, 483, 7, 141, 246 and international counterparts. Other patents pending.

M-UR-Vaqua/GB-07/2007

 **COVIDIEN**

COVIDIEN UK LTD
154 FAREHAM ROAD
GOSPORT HAMPSHIRE PO13 0AS
UNITED KINGDOM

WWW.COVIDIEN.COM
WWW.VAQUACATHETER.COM

**New VaQua coated intermittent catheter.
Makes life easier than ever for your patients.**

VaQua is a new, coated intermittent catheter from Covidien, a company with wide urological experience. Drawing on our knowledge in this area, we have developed and designed VaQua coated intermittent catheter around

the comfort of the people who use it. As such, it incorporates a number of medically and technologically advanced features that make the self-catheterisation process easier and more comfortable than it has ever been.

5-SECOND ACTIVATION

VaQua coated intermittent catheter has a catheter activation time of just 5 seconds! In contrast, almost all other hydrophilic intermittent catheters must be immersed in water for at least 30 seconds to activate the coating. In this way, VaQua coated intermittent catheter will help your patients to get on with their daily lives, however busy they might be.

DISCREET

The packaging is discreet, waterproof and easy to open.

SMOOTHER DRAINAGE EYES

The drainage eyes are made by an ultrasonic method. This ensures maximum smoothness and added comfort.

PVC CATHETER MATERIAL

VaQua coated intermittent catheter is constructed from a high quality medical grade PVC. It is rigid enough to aid insertion, yet flexible enough to prevent trauma and discomfort.

SMOOTH, ROUNDED CATHETER TIP

Ease and comfort of insertion is further enhanced by a catheter tip that is designed to reduce the risk of trauma.

SMOOTHER COATING FOR MINIMAL FRICTION

The patented hydrophilic coating binds with the activation fluid to create a surface that is essentially frictionless. This has two benefits during catheter insertion and removal:

- the catheterisation process becomes more comfortable,
- the risk of urethral trauma is significantly reduced.

The right catheter for your patients

VaQua coated intermittent catheter is available in Male, Female and Child versions:

- Male: 40 cm
- Female: 15 cm
- Child: 30 cm

Within those versions, VaQua coated intermittent catheter is available in a wide range of diameters from 8 Fr/Ch to 18 Fr/Ch.

By enabling you to select from a wide range of sizes, VaQua coated intermittent catheter offers the highest standards of comfort and convenient self-catheterisation for all your patients. Colour coding of the catheter funnel allows you to quickly identify the catheter size. If you would like a free sample kit so that your patients can try VaQua coated intermittent catheter, please contact your local Covidien representative.



Product ordering information

	Code	Description	Pack Size
PAEDIATRIC 30 CM	VAQ08P	8 Fr/Ch	30
	VAQ10P	10 Fr/Ch	30
FEMALE 15 CM	VAQ08F	8 Fr/Ch	30
	VAQ10F	10 Fr/Ch	30
	VAQ12F	12 Fr/Ch	30
	VAQ14F	14 Fr/Ch	30
	VAQ16F	16 Fr/Ch	30
	VAQ18F	18 Fr/Ch	30
MALE 40 CM	VAQ08M	8 Fr/Ch	30
	VAQ10M	10 Fr/Ch	30
	VAQ12M	12 Fr/Ch	30
	VAQ14M	14 Fr/Ch	30
	VAQ16M	16 Fr/Ch	30
	VAQ18M	18 Fr/Ch	30